



arplay medical



Storage furniture for treatment room in radiotherapy

By taking your comments and requirements into account, ARPLAY has developed a customized solution which can be adapted to any treatment room you have available.

This furniture is designed to store your electron cones with their shieldings, your lead trays, breast boards, masks, head rests and all the accessories you have.

The furniture is made from laminate panels, 19mm thick, with PVC edgings.

The drawers with PVC edges, are mounted on 50kg load bearing runners and, as there are buffers at the end of these runners, these drawers can be pulled right out.

The doors have PVC edges and are mounted on invisible 110° opening hinges. The furniture is mounted on 4 jacks to level it. The choice of colours is made from a colour chart that we provide in order to obtain a perfect integration of the furniture in its environment.



Shelves for lead trays and foam blocks

Shelves for masks and shieldings

